



[www.hearingtoday.org](http://www.hearingtoday.org)

## A Quick Reference Guide to Hearing Aids

### For mild to moderate hearing loss

**Invisible In The Canal** is a custom-fitted style inserted into the ear canal making it completely invisible when worn.



**Completely in Canal** are custom-made to fit in the ear canal. Only the tip of a small handle shows which is used to insert and remove the instrument.



### For mild to mildly severe hearing loss

**In the Canal** hearing instruments are custom-made to fit in the ear canal with a smaller portion of the device showing in the outer ear.



### For mild to severe hearing loss

**In the Ear** hearing instruments are custom-made to fit within the outer portion of the ear.



**Receiver in Canal** has the receiver inside the ear canal. Thin electrical wires are used to reduce distortion.



### For mild to profound hearing loss

**Behind the Ear** devices are housed in a casing that rests behind the ear. A clear plastic tube directs sound into an earmold fitted inside the ear canal.

